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Dear customer,

Before operating your new appliance, please read these operating instructions carefully. They contain important information for safe use, for installation and for maintenance of the appliance. Please keep these operating instructions for future reference. Pass them on to possible new owners of the appliance.

Notes: which are important for your safety or for the proper functioning of the appliance, are stressed with a warning triangle and or with signal words (**Warning! Caution! Important!**). Please observe the following carefully.

1. GENERAL SAFETY INSTRUCTIONS

For your own safety and for the proper utilization of the appliance, read this manual carefully, including the warnings and recommendations, before installing the appliance and using it for the first time. In order to avoid damage to the appliance and/or personal injury, people using the appliance should be fully familiar with its operation and safety functions. Keep these recommendations close to the appliance, for future reference, and ensure that the document is transferred with the appliance, if it is sold or in the case of moving. This will ensure optimal performance of the appliance.

To avoid any risk of injury, keep this manual. The manufacturer declines any liability whatsoever in the case of operating/handling errors concerning the appliance.

Safety of children and other vulnerable people

- This appliance may only be used by specialists and laboratory personnel who have been trained for this purpose and are familiar with all safety measures pertaining to laboratory work. **Children and people with reduced physical, sensory or mental capabilities or lack of experience and knowledge must not commission or operate this appliance.**
- Keep all packaging out of the reach of children, as there is a risk of suffocation.
- If you decide to scrap the appliance, disconnect it from the wall socket, cut the connection cable (as close to the appliance as possible) and remove the door, in order to prevent children from playing with it, with the resulting risk of electric shock or being trapped inside.
- If an appliance with a magnetic door seal is being replaced by another appliance with a latch on the door or on the seal, ensure that the latch is disabled before selling or scrapping the old appliance. This will prevent the appliance from being transformed into a trap for children.
- Keep the appliance and its power cable out of reach of children less than 8 years old.
- Medical appliances can only be operated by trained and authorized personnel

General safety

CAUTION — Never store inflammable substances, such as aerosols, inside the appliance, as they could leak or release fluids.

CAUTION — A damaged power cable should be immediately replaced by the manufacturer or by a qualified retailer, in order to avoid any risk of electric shock.

CAUTION — Keep the ventilation openings well open and ensure that air can circulate.

CAUTION — Do not damage the refrigerating circuit.

CAUTION — Do not use electrical appliances inside the compartment, unless they are tolerated by the manufacturer.

Caution: Fire risk

If the refrigerating circuit is damaged:

- Avoid open flames or any source of combustion.
- Ensure that the room in which the appliance is installed is well ventilated.
- Alcoholic liquids or other containers with alcohol inside must only be stored tightly sealed. Any alcohol that leaks out may be ignited by electrical components.

It is dangerous to change the specification of the appliance in any way whatsoever. Damage to the power cable may cause a short circuit and/or electric shock.

Electrical safety

- We decline any liability in respect of incidents caused by poor electrical installation.
- Do not extend the power cable. Do not use an extension cable, adapter, or multi-way socket.
- Check that the wall socket is not damaged. A wall socket in poor condition may cause the appliance to overheat and explode.
- Check that the wall socket used for the appliance is properly accessible and may be unplugged in case of emergency.
- Never pull on the main power cable.
- Ensure that the appliance is always electrically earthed.
- If the wall socket is loose, do not connect the cable, as there is a risk of electrical shock or fire.
- Do not use the appliance if the interior lighting cover is not installed.
- Disconnect the appliance before replacing the LED lighting which should be replaced by a professional.
- The appliance operates on a single phase 220~240 V / 50 Hz power supply. The appliance must be connected to an earthed wall socket, in accordance with applicable recommendations.
- Do not attempt to replace the damaged power cable, contact your after-sales service.
- The power socket must be easily accessible but out of the reach of children. In case of doubt, contact your installation technician.

Range of use of appliance

Correct use

- Intended use This medical refrigerator is suitable for professional storage of products at a temperature of 5°C.
- Typical products to be stored include drugs which require refrigeration, medicines or other pharmaceutical products in a temperature range from 2°C to 8°C.
- The appliance meets the requirements of DIN 13277 Refrigerators and freezers for laboratory and medical applications.

- It is necessary to use an independent alarm system with permanent monitoring if storing valuable or temperature-sensitive substances or products. This alarm system must be designed in such a way that every alarm state is immediately registered by a competent person who then takes appropriate measures.

Foreseeable incorrect use

Do not use the appliance for the following applications:

- Storage and refrigeration of:
 - Chemically unstable, flammable or corrosive substances
 - Blood, plasma or other body fluids for the purpose of infusion, application or introduction in the human body.
- Use in potentially explosive atmospheres.
- Use outdoors or in areas exposed to dampness and splashing water.

Improper use of the appliances will result in goods stored being damaged or spoiled.

Daily utilization

The appliance may only be installed, tested, maintained, and commissioned by specialist personnel who are familiar with the installation, commissioning, and operation of the appliance. Specialist personnel are people who, on account of their specialist training, knowledge and experience as well as their knowledge of the relevant standards, are able to assess and perform the work assigned to them and identify potential hazards. They must have training, instruction, and authorization to work on the appliance

- Do not store inflammable substances or fluids inside the appliance. This would create a risk of an explosion.
- Do not operate other electrical equipment inside the appliance.
- When disconnecting the appliance, grip the plug firmly and do not pull on the cable.
- Do not expose the appliance to direct sunlight.
- Keep the appliance away from candles, luminaires and open flames, in order to avoid fire risk.
- This appliance is only designed to store medicine.
- The appliance is heavy. Be careful when moving it.
- If your appliance is equipped with castors, remember that these are only designed to facilitate minor movements. Do not move the appliance over long distances.
- Do not lean or climb on the appliance or its elements.
- In order to avoid falling objects and to prevent damage to the appliance, do not overload the compartments of the appliance.

Danger of frostbite, numbness and pain:

Avoid prolonged skin contact with cold surfaces or refrigerated/frozen goods or take protective measures, e.g. wear gloves.

Risk of crushing

Do not hold the hinge when opening and closing the door. Fingers may get caught.

Caution! Cleaning and servicing

- Before carrying out servicing, switch off the power supply and disconnect the appliance.
- When cleaning the appliance, do not use metal objects, steam systems, volatile fluids, organic solvents or abrasive substances.
- Do not use sharp or pointed objects to remove ice. Use a plastic scraper.

Important information concerning installation!

- To ensure good electrical connection, follow the recommendations given in the present manual.
- Unpack the appliance and visually check that it is not damaged. Do not connect to a damaged appliance. Inform your retailer or point of sale of any damage. In such a case, keep the packaging.
- It is recommended that the appliance be allowed to settle for at least 4 hours before connecting the electrical power supply, so that the compressor fluid feed is correct.
- Clean air circulation is necessary to avoid overheating. To ensure sufficient ventilation, observe the installation recommendations provided.
- The more refrigerant there is in the appliance, the larger the space it is installed in must be. If the space is too small, any leak may create a flammable mixture of gas and air. For every 8 g of refrigerant, the installation space must be at least 1 m³. Information regarding the refrigerant can be found on the type of plate inside the appliance.
- The floor of the location must be horizontal and even.
- The installation location must be able to withstand the weight of the appliance plus weight when stocked to maximum level.
- In order to avoid fire risk, wherever possible, ensure that the appliance is not in contact with walls or with any hot element (compressor, condenser). Observe the installation recommendations.
- The appliance should not be placed close to radiators or cooking hot plates or gas rings.
- Check that electrical sockets remain accessible when the appliance is installed.
- During transportation or use, no violent vibration or collision is allowed, and the appliance shall be kept away from rain. Store in a clean room with humidity no more than 80 %, no corrosive gas and good ventilation

POWER SAVING

To limit the electricity consumption of your appliance:

- Install it in an appropriate place (see “Installing your appliance”).
- Leave the doors open as little time as possible.
- To ensure that the appliance operates properly, clean the condenser regularly (see “CARE AND MAINTENANCE” chapter).
- Regularly inspect the door seals and ensure that they are still closing properly. If this is not the case, contact your after-sales service.

Troubleshooting

- Any electrical work should be carried out by a qualified and skilled technician.
- The appliance should be repaired by an accredited repair center, using only manufacturer's original spare parts.

The spare parts list is available from the manufacturer; we invite you to contact MEDIFROID for any need of spare parts:

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R600A

Safety instructions

Caution: Do not obstruct the ventilation grills on the appliance.

Caution: Do not store explosive substances, such as aerosol cans containing an inflammable propellant, in the appliance.

Caution: Do not damage the refrigerating circuit in the appliance.


Caution: The refrigerating units use isobutane (R600a) and should not be placed close to ignition sources (for example, electrical contacts that are exposed or open contacts that could be closed by the refrigerating fluid in the event of a leak). The refrigerating fluid type is specified on the power cabinet identification plate.

Caution: Do not use electrical equipment in the compartments of the appliance, unless it is of the type recommended by the manufacturer.

This appliance meets the requirements of all applicable European directives and their amendments.

R290

Safety instructions

Caution: Do not obstruct the ventilation grills on the appliance.

Caution: Do not store explosive substances, such as aerosol cans containing an inflammable propellant, in the appliance.

Caution: Do not damage the refrigerating circuit in the appliance.


Caution: The refrigerating units use isobutane (R290) and should not be placed close to ignition sources (for example, electrical contacts that are exposed or open contacts that could be closed by the refrigerating fluid in the event of a leak). The refrigerating fluid type is specified on the power cabinet identification plate.


Caution: Do not use electrical equipment in the compartments of the appliance, unless it is of the type recommended by the manufacturer.

This appliance meets the requirements of all applicable European directives and their amendments.

2. PRODUCT DESCRIPTION AND DIMENSIONS

PHAR200


PHAR300 & PHAR400

Model	Loading capacity (shelves)	Loading capacity (drawer)
PHAR200	21kg	10kg
PHAR300	20kg	10kg
PHAR400	20kg	10kg

Dimensions in mm.

		PHAR200	PHAR300	PHAR400
A	Height	810	1762	1992
B	Length	650	650	650
C	Depth	625	673	673
1	Depth	625	673	673
2	Length	650	650	650
3	Length with open door	1092	1098	1098
4	Depth with open door	1060	1084	1084

3. INSTALLATION INSTRUCTIONS

Location

When selecting a position for your appliance, you should make sure the floor is flat and firm, and the room is well ventilated. Avoid locating your appliance near a heat source, eg, cooker, boiler or radiator. Also, avoid direct sunlight in out-buildings or sun lounges. Avoid any very humid environment, otherwise condensation will occur on the outside of the appliance cabinet. Never place the appliance on a recessed wall or into fitted cabinets or furniture when your appliance is working, the grill at the back may become hot and the sides warm. **It must therefore be installed so that the back and the sides of the appliance have at least 30 cm of free space.**

In order to improve efficiency of the cooling system and save energy, it is necessary to maintain good ventilation around the appliance for the dissipation of heat. For this reason, sufficient clear space should be available around the appliance.

Climate class

When selecting a position for your appliance, you should make sure the floor is flat and firm, and the room is well ventilated. An average ambient temperature must be observed.

This temperature is specific according to the climatic class of the device:


CLASS	SYMBOL	AMBIENT TEMPERATURES RANGE (°C)
Extended temperate	SN	From +10 to +32
Temperate	N	From +16 to +32
Subtropical	ST	From +16 to +38
Tropical	T	From +16 to +43

N: This refrigerating appliance is intended to be used at ambient temperatures ranging from 16 °C to 32 °C

Installing multiple appliances

Risk of damage due to condensation between the side walls.

- Do not install the appliance directly next to another fridge.
- Install appliances with a space 10 cm between appliances.
- At higher levels of humidity, increase space between


Preparation before use:

1. Remove the outer packaging of all products (including the protective foam in the packaging box)

Warning: Don't put the plastic bags within the reach of children, so as to prevent suffocation accidents.

2. Inventory of accessories: Please check the accessories and materials according to the packing list.
3. Cleaning: Clean the product once before use (check the information below).
4. Cabinet fixing: twirl the support leg or lock the universal caster to ensure the cabinet is not moving when in use.
5. Install the fixing plates: this appliance is supplied with a safety spacer which allows sufficient distance to be maintained between the appliance and the wall:
 - A. Firstly, take out the two fixing plates and check if they are consistent
 - B. Remove the bolts from the fixing plate and then fix the plate on the back of the refrigerator.
 - C. Push the appliance with the mounting parts mounted against the wall.

Warning : If you accessorize your appliance with drawers, be careful to only open one drawer at a time to avoid your appliance tipping over.

Leveling the appliance

If the appliance is not leveled, the door and magnetic seal alignment will be affected and may cause your appliance to work incorrectly. Once the appliance is placed in its final location, adjust the leveling feet at the front by turning them.

Cleaning before use

The equipment should be cleaned and maintained regularly (for safety, please unplug the power plug), and the inner and outer surfaces of the cabinet should be wiped with a warm, damp soft cloth.

Anti-tip system (PHAR200 - PHAR300 – PHAR400)

The PHAR300 and PHAR400 showcases are equipped with an anti-tip system. To install the anti-tip system, you will need a screwdriver and an electric drill.

1. Screw the 2 metal brackets into the holes provided on the back of the showcase. Secure each bracket with 4 round-head screws.

2. Mark the location of the metal brackets (height and spacing) on the wall.

3. Drill the mounting holes in the wall at the previously marked locations.

4. Use a hammer to place the provided expansion screws into the drilled holes.

5. Position the showcase at the expansion screws and slide the unit to the right so that the screws are positioned at the two previously installed metal brackets. The screws hold the unit in place and prevent it from tipping over.

6. Use a suitable wrench to tighten the expansion screw bolts and secure the display case. Once the bolts are tightened, the display case is securely held in place.


4. OPERATING YOUR APPLIANCE

First Power-on


When using the equipment for the first time, please follow these steps:

1. After the equipment is placed, leveled and cleaned, it shall stand for more than 24 hours, and then power it on to ensure the normal operation of the equipment.
2. Under no-load condition, connect the power cord to a special outlet with appropriate specifications.
3. **After powering on, turn on the power switch of the equipment. Toggle the on/off battery switch to the ON position, located in the back of the unit.**

PHAR200

PHAR300 & PHAR400


4. When you first power-on the appliance, the alarm (high temperature alarm) will be activated (as interior appliance temperature is higher than the default adjusted one 5°C). Press on mute button  to stop the buzzer each time it is activated till the interior appliance temperature will match the set value.

5. The temperature inside the appliance is set to 5°C (+/- 3°C) as required DIN 13277 and we recommend not to change it.

Note: Please contact us for any temperature setting requirements.

6. Set the date and the hour (time settings) as following :

- Use the " or "
" key to navigate through parameters till you reach the parameter "n" (to set the year).
- Press the "
" key to display the parameter value.
- Use " or "
" keys to adjust the parameter values.
- Press the "
" key to temporary store changes and return to the parameter list
- Repeat the process for other date and hour parameters :

Setting	Description	Value
n	Set logger module time - year	10~50
y	Set logger module time - month	1~12
r	Set logger module time - day	01~31
s	Set logger module time - hour	00~23
F	Set logger module time - minute	00~59

- f. To save all modifications, press and hold the "✓" key for more than 3 seconds.
- g. To exit the parameter setting mode, either press the "✓" key for more than 3 seconds or **wait 60 seconds without pressing any key.**

Note: If the digital display flashes "LoF", the recorder is not activated. To activate, press and hold both ✓ and ▲ keys for 3 seconds until "LoF" disappears.

WARNING: Please do not change any other setting, for any other requirements, contact your after sales service.

Storage

1. Check whether the operating temperature of the equipment reaches the required value, observe the normal start and stop of the equipment for more than 24 hours, and put a small number of items in the appliance after confirming the normal performance.
2. Please store items in batches, with the items not exceeding 1/3 of the cabinet volume each time. Ensure that the equipment is running properly after a shutdown for more than 12 hours before putting in the next batch of items).
3. Observe maximum load, find below the maximum loading capacity for each

model . 

Model	Loading capacity (shelves)	Loading capacity (drawer)
PHAR200	21kg	10kg
PHAR300	20kg	10kg
PHAR400	20kg	10kg


WARNING: Do not stack refrigerated products above the indicated height. This is important for the air circulation and even temperature distribution in the interior of the product.

4. Keep liquids in closed containers
5. Try not to open the door during cooling, otherwise the temperature will rise.

Control panel:

This series of products includes two control systems, based on which different adjustment methods apply. Please select the proper adjustment method according to the type of control system with the product that you purchased.

PHAR200

a) Temperature Display Window

Displays the average internal temperature of the cabinet in °C under normal operation.

b) Door Open Indicator

- Lights up when the door is open.
- If the door remains open for more than 1 minute, an alarm is triggered, and "do" is displayed.

c) Alarm Indicator

- Stays off when the product operates normally.
- Lights up when an operational abnormality occurs.

d) Light Control

- After powering on, the light is off by default.
- The light can be turned on or off using the designated key.

e) Down Key :

- In parameter setting mode, parameter values decrease (e.g., lowering the set temperature).
- Holding the button speeds up the decrease.

f) Up Key :

- In parameter setting mode, increases parameter values (e.g., raising the set temperature).
- Holding the button speeds up the increase.
- Holding for 3 seconds under normal operation imports USB data from the past 12 months.

g) Set/Mute Key :

- In case of no alarm state and key unlock state, press , and display the ambient temperature for 5s and then return to normal display.
- In the unlock state, press  for more than 3s, and enter the user setting menu.
- When the buzzer is triggered (including cabinet high temperature alarm, door opening alarm, sensor failure alarm, etc.) and in the key unlock state, press  for the first time, and the buzzer stops ringing, and the ambient temperature is displayed for 5s, after which the normal display is resumed (pressing the mute button is only to turn off the buzzer for alarming this abnormal state, for example troubleshooting, and the buzzer will be triggered next time for any abnormality). Then press  again, trigger the buzzer, display the ambient temperature for 5s.

h) USB Interface

Automatic USB Data Export

- When a USB drive is connected, the recorder beeps once and displays "ON".
- The system generates PDF files of the data from the current and previous months.
- After completion, the recorder beeps again, displays "End", and returns to normal display after 6 seconds.

Note: If data is minimal, "ON" and "End" may not appear.


Manual USB Data Export

- In the unlocked state, if a USB drive is connected and no file is currently being generated:


- Press and hold the “▲” key for 3 seconds; the display will show "d01".
- Use the “▲” or “▼” keys to adjust between "d00" to "d12".
- Press the designated key to cancel file generation (d00) or create a PDF file of data from previous months (1-12).

Note: If the digital display flashes "LoF", the recorder is not activated. To activate, press and hold both  and  keys for 3 seconds until "LoF" disappears.

PHAR300 & PHAR400

1. Indicators description

Door switch	WIFI	Defrosting	Key lock
Power failure	Print (option)	Door heating	Refrigeration
Mute	Electric quantity	Serial port	Fan

a) Door Open Indicator

- Lights up when the door is open.
- Turns off when the door is closed.

b) WIFI Indicator

This function is not available on this model.

c) Defrost Indicator

- Remains on when the display case is in defrost mode.
- Turns off once defrosting is complete.

d) Key Lock Indicator

- When the keypad is locked, no keys respond, and the indicator lights up.
- To unlock, press and hold the "Up" + "Down" keys for 3 seconds. You will be prompted to enter the default password "005." Confirm with the "Set/Mute Multiplex" key to disable the lock (the indicator turns off).
- In the unlocked state, if no key is pressed for 60 seconds, the lock activates automatically, and the indicator lights up.
- To manually lock, press and hold the "Up" + "Down" keys for 3 seconds.

e) Power Failure Indicator

- The device normally operates on 220V/110V.
- In case of a power outage, a buzzer is activated, the display alternates flashing the code "PF" every 3 seconds, and the power failure indicator lights up.
- When the power is restored, normal operation resumes, and the indicator turns off.

f) Print Indicator

This function is not available on this model.

g) Door Heating Indicator

- Lights up when door heating is activated.
- Turns off when door heating is deactivated.

h) Refrigeration Indicator

- Lights up when the compressor is running.
- Turns off when the compressor is stopped.

i) Mute Indicator

- Lights up when the alarm sound is muted.
- Turns off when the alarm sound is reactivated.

j) Low Battery Indicator

- When the battery voltage is below 8V, a buzzer sounds, the low battery indicator lights up, and the display alternates flashing the code "BL" every 3 seconds.
- When the battery voltage exceeds 12V, the buzzer stops, the indicator turns off, and the display returns to normal.

k) Serial Port Indicator

- Turns off when the RS-485 serial port is not connected.
- Lights up when the equipment is correctly connected to the RS-485 serial port.

l) Fan Indicator

- Lights up when the evaporator fan is running.
- Turns off when the evaporator fan is stopped.

2. Temperature display window:

Which displays the average temperature inside the cabinet in °C under normal operation.

Temperature check:

- In the key lock state, press  key, and the control panel displays the ambient temperature, and returns to normal display after 5 s without key operation, or pressing  and .
- In the key unlock state, press  key, and the control panel displays ambient temperature, and returns to normal display after 5s without key operation.

Humidity check:

key unlocked state, long press  and  , digital tube display humidity, no press any Key operation after 5 seconds or press  and , return to normal display.

3. Set/Mute key :

- In case of no alarm state and key unlock state, press , and display the ambient temperature for 5s and then return to normal display.
- In the unlock state, press  for more than 3s, and enter the user menu.
- When the buzzer is triggered (including cabinet high temperature alarm, door opening alarm, sensor failure alarm, etc.) and in the key unlock state, press  for the first time, the buzzer stops ringing, and the ambient temperature is displayed for 5s, after which the normal display is resumed (pressing the mute button is only to turn off the buzzer

for alarming this abnormal state, for example troubleshooting, and the buzzer will be triggered next time for any abnormality). Then press  again, trigger the buzzer, display the ambient temperature for 5s.

- In unlock state and parameter setting mode, press this key to display parameter values and parameter names. If the pressing time is longer than 3 seconds, it saves the settings and returns to the normal interface.

4. **Up key**

- In parameter setting mode, move to the next parameter or increase the parameter value. For example, when setting the set temperature increases the set temperature value.
- When setting the parameter value, long press the “Up” button, and the parameter will increase rapidly.
- Under normal conditions, long press the “Up” key for 3 seconds to import the data of the USB flash drive in 12 months.

5. **Down key**

- In parameter setting mode, move to the previous parameter or decrease the parameter value. For example, when setting the set temperature, reduce the set temperature.
- When setting the parameter value, long press the down button, and the parameter will decrease rapidly.

6. **Light switch key**

When the equipment is powered on, the lights are turned off by default.

7. **USB data export:**

- Automatic USB Data Export
 - When a USB drive is connected, the recorder beeps once and displays "ON".
 - The system generates PDF files of the data from the current and previous months.
 - After completion, the recorder beeps again, displays "End," and returns to normal display after 6 seconds.

Note: If data is minimal, "ON" and "End" may not appear.

- Manual USB Data Export
 - In the unlocked state, if a USB drive is connected and no file is currently being generated:
 - Press and hold the “Up” key for 3 seconds; the display will show "d01".
 - Use the “Up” or “Down” keys to adjust between "d00" to "d12".
 - Press the designated key to cancel file generation (d00) or create a PDF file of data from previous months (1-12).

Note: If the digital display flashes "LoF," the recorder is not activated. To activate, press and hold both  and  keys for 3 seconds until "LoF" disappears.

Function Setting (all models):

- After powering on, the equipment can enter the working state.
- User parameter settings:
 1. **Unlocking:** In normal operation, press and hold the  +  keys for 3 seconds. The display will show "000".
 2. Enter the password "005" using the  key and confirm by holding the  key to unlock.

3. Press the "✓" key for 3 seconds. The display will show the different parameters to be adjusted.
4. Use the "▲" or "▼" key to navigate through parameters.
5. Press the "✓" key to display the parameter value.
6. Use "▲" or "▼" keys to adjust the parameter values.
7. Press the "✓" key to temporary store changes and return to the parameter list
8. Repeat the process for other parameters as needed.
9. To save all modifications, press and hold the "✓" key for more than 3 seconds.
10. To exit the parameter setting mode, either press the "✓" key for more than 3 seconds or **wait 60 seconds without pressing any key.**

Parameter display

Menu item	Parameter Range	Suggested settings	Remarks
MAX	—	—	The highest temperature since last clearance
MIN	—	—	The lowest temperature since last clearance
CLR	—	—	Clearance of the Max and Min temperature records
n	Set logger module time - year	10~50	
y	Set logger module time - month	1~12	--
r	Set logger module time - day	01~31	--
s	Set logger module time - hour	00~23	--
F	Set logger module time - minute	00~59	--

If there is no operation for 60 seconds under the quick setting menu, it will automatically exit the quick setting menu and return to normal display.

Alarm display error codes (all models):

Code	Error Description
H1	High temperature alarm
L1	Low temperature alarm
H2	Alarm for high ambient temperature
H3	Condenser overheat alarm
do	Door opening alarm
PF	Power failure alarm
bL	Battery low alarm
Er	The recorder is not connected
LoF	Recorder did not start
EE	Communication failure

Note: This appliance is equipped with a “**safety thermostat function**”; in case that interior appliance temperature goes below 0°C, the compressor will stop automatically.

Alarm demonstration function (all models):

Under "key lock state" Long press  the digital tube displays the temperature value and alarm code of the over-temperature alarm, and the buzzer beeps.

Maximum and minimum temperature display (all models):

Under the key lock state:

- Short press the  to view the maximum value: after blinking 3 times, the maximum value menu and the maximum temperature displays.
- Short press  to view the minimum value, after blinking 3 times, the minimum value menu and maximum temperature displays.
- Long press combination key 3 seconds  +  +  to clear the current maximum and minimum value data, at the same time, the buzzer beeps a reminder.

Notes:

1. **Memory function:** The appliance retains the set values. After a power failure, it resumes operation with the last configured settings.
2. **Compressor protection:** If the appliance is turned off, wait at least 5 minutes before turning it back on to prevent compressor damage.
3. **Storage responsibility:** The equipment ensures normal operation under specified conditions, but the manufacturer is not liable for loss or damage to stored items during power failure.
4. **Monitoring and emergency handling:** A designated individual should check and record the appliance's operational status every 2–4 hours. In case of failure or shutdown, the internal temperature will rise. If repairs are delayed, transfer stored items to another facility that meets the required storage conditions.
5. **Pre-storage temperature check:** Verify that the appliance's temperature range meets the storage requirements of the items. Ensure that neither the air inlet nor outlet is obstructed when placing items inside.
6. **Temperature fluctuations:** Due to appliance inertia, slight differences between the displayed temperature and the set temperature are normal. Avoid frequent changes to the set temperature within a short time, as large temperature fluctuations may reduce efficiency.
7. **Item placement guidelines:** Items should not be placed directly on the cabinet's bottom; always use a bottom shelf for optimal refrigeration performance.
8. **Humidity Considerations:**
 - Excessively moist or dry items can impact internal humidity levels.
 - To maintain optimal conditions, store items in sealed containers.
 - Frequent or improper door usage affects internal humidity stability.


Following these guidelines ensures optimal performance and longevity of the appliance while maintaining the safety and integrity of stored items

Other equipments:

▪ Remote alarm terminal and RS485 interface (all models):

Your device is equipped with a remote alarm connection to send alarm signals to a building management system, a monitoring unit, visual or acoustic warning device etc. The alarm output is managed by the built-in controller of the unit.

Alarm criteria may be too low temperature, too high temperature, power failure, open door or similar.


The alarm output is an electrical relay which acts like an ON/OFF switch.

▪ Hole test

25 mm Portholes are used to pass the temperature probe inside the cavity


1. Remove the two rubber plugs.
2. Press a pencil through the marks on the inside of the rubber plugs.
3. Feed your external sensor through the passage.
4. Feed the sensor through the outer rubber stopper

▪ Test bottle

The ratio of the mixed liquid is 90ml of pure water and 10ml of glycerol.

The mixture in the probe bottle shall be checked and refilled on a regular basis.


▪ Sensors


Appliances are equipped with the following sensors for temperature monitoring

PHAR300 & PHAR400

Seven Sensors, Precise temperature control—Multiple probe distribution

PHAR200

5. CARE AND MAINTENANCE

Defrost

Your appliance is equipped with an automatic defrost cycle. When a cooling cycle finishes, the appliance's refrigerated surfaces are defrosted automatically. The defrost water is channeled into a condensates' evaporation tray which is located in the back of the appliance near the compressor. The heat produced by the compressor then evaporates the condensates collected in the tray.

Cleaning inside the appliance

Before maintenance, switch off the appliance and disconnect the mains plug from the mains socket. After defrosting, you should clean the appliance internally with a weak solution of bicarbonate of soda. Then rinse with warm water using a damp sponge or cloth and wipe dry. Wash the baskets in warm soapy water and ensure they are completely dry before replacing

them in the appliance. Condensation will form on the back wall of the appliance. However, it will normally run down the back wall and into the drain hole behind the salad bin.

Cleaning outside the appliance

Before maintenance, switch off the appliance and disconnect the mains plug from the mains socket. Use standard non-abrasive detergent diluted in warm water to clean the appliance exterior.

The grill of the condenser at the back of the appliance and the adjacent components can be vacuumed using a soft brush attachment.

Do not use harsh cleaners, scouring pads or solvents to clean any part of the medical showcase.

Condenser cleaning: clean the dust on the condenser regularly to avoid affecting the cooling effect (cleaning method: first open the screws of the ventilation hood or front hood and use a brush to gently brush away the dust on the surface of the condenser).

Moving the unit

▪ **Location**

Do not place your appliance near a heat source, for example. Stove, boiler or radiator. Also, avoid direct sunlight.

▪ **Showcase upgrade**

Make sure that the showcase is level. Use the rotatable adjustable feet at the front. If the showcase is not level, the doors and magnetic sealing alignments will be affected and may cause malfunctions.

▪ **Do not turn on the showcase for 4 hours**


After any movement, it must not be connected for 4 hours. This delay allows the coolant to be adjusted.

▪ **Installation**

Let your device breathe. Heat is released from the walls, so it is important to leave a circulation of air around your device.

Changing the LED light

This appliance consists of light-emitting diodes (LEDs). This type of diode cannot be changed by the consumer. The lifetime of these LEDs is normally sufficient for no change to be made. If however, and despite all the care taken during the manufacture of your appliance, the LEDs were defective, please contact your after-sales service for any intervention.

6. RECHARGEABLE BATTERIE

This appliance is equipped with a battery backup system, which supplies the controller with power at power failure. The backup system duration is 48 hours.

The battery is located in the top of the appliance, inside the electrical box.

Battery maintenance

- a. To prolong battery life and prevent prolonged inactivity, operate the product for at least 24 hours each month to ensure proper charging.

- b. If the appliance is unused for an extended period, connect it to a power supply at least once a month. Turn on both the power switch and battery switch to complete the charging process. Each charging session should last at least 24 hours.
- c. After discontinuing equipment use, turn off both the power switch and battery switch promptly. Failure to do so may result in battery capacity loss and, in severe cases, permanent battery damage.
- d. During long-term power outages or transportation, always turn off the main power switch. Failure to do so may lead to prolonged battery discharge, power loss, and even permanent battery damage, potentially causing display malfunctions upon reactivation.
- e. The battery is a consumable component with an expected lifespan of approximately 2-3 years. Improper usage, such as allowing power loss or reaching the end of its service life, may trigger a low battery alarm. While refrigeration remains unaffected, the alarm will persist.

Battery replacement and recovery

The electrical box contains high-voltage components. Only qualified engineers or maintenance personnel should open the electrical box. So, please contact after-sales support to change the battery.

Note: For a replaced battery that can be recycled, please contact the local battery recycling agency.

7. TROUBLESHOOTING GUIDE

Any product may fail. We recommend yearly maintenance by a professional.

Please observe the operation of equipment in time during use. If there is any abnormality, please check and handle it according to the following table first. If the abnormality can't be changed, please inform our service center in time, and we will serve you wholeheartedly to avoid losses.

Other problems	Causes and solutions
Equipment does not work	Please make sure that the outlet is energized. Please make sure that the power plug is plugged in, not loose. Please make sure the power fuse is not disconnected. Please make sure the supply voltage is appropriate, not too low or too high.
Compressor is not running	Please make sure that the temperature is set correctly. Please check whether the temperature inside the cabinet is too low.
The temperature does not reach the set value	Please make sure that the door is closed tightly and don't open it too many times during a short time. Please don't put too many items in at one time. Please make sure that the ambient temperature is not too high.
High noise	Please make sure that the cabinet is placed on a flat ground. Please make sure that the cabinet does not contact the wall.
Condensation on cabinet surface	In rainy and humid seasons, door condensation is normal, and it shall be wiped off with a dry cloth.
Undesirable odor	Equipment needs to be cleaned. Items with heavy smell are not packed.
The door is not closed properly, and the cool air leaks	Please make sure that the temperature is set correctly. Please check whether the temperature inside the cabinet is too low.

The alarm lamp flashes, and the buzzer is triggered	After the equipment is used for a period of time, the door seal becomes hard and deformed. Maintenance method: Blow the deformed part of the hot seal with a blower to soften it, and then close and compress it after the door seal becomes soft.
Lighting damage	Please contact after-sales service personnel for the replacement. Please do not replace the parts by yourself.

Alarm display error codes: contact after sales service personnel

Code	Error Description
H1	High temperature alarm
L1	Low temperature alarm
H2	Alarm for high ambient temperature
H3	Condenser overheat alarm
do	Door opening alarm
PF	Power failure alarm
bL	Battery low alarm
Er	The recorder is not connected
LoF	Recorder did not start
EE	Communication failure

The following conditions are not faults

Compressor Noise

When the compressor starts or stops, slight impact sounds from the equipment parts may occur. This is a normal noise.

High Temperature and High Humidity Alarm

If overheated items are placed inside after opening the door, the control system may trigger high temperature and high humidity alarms (if this function is available, refer to the alarm display table).

Solutions:

- Allow items to cool naturally before placing them inside.
- Load items in small quantities and in batches rather than all at once.
- The alarms will automatically deactivate once the system stabilizes.

Refrigerant Flow Sound

A slight sound similar to running water may be heard due to refrigerant flowing through the pipeline, which is normal during operation.

8. ENVIRONMENT

This appliance is marked according to the European directive 2012/19/CE on Waste Electrical and Electronic Equipment (WEEE). By ensuring this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste handling of this product.

The symbol on the product, or on the documents accompanying the product, the Crossed Out Wheelie Bin, indicates that this appliance may not be treated as household waste. Instead, it shall be handed over to the applicable collection point for the recycling of electrical and electronic equipment, or, if you are buying an equivalent product, to the retailer of the new product.

The user is responsible for taking the appliance to the appropriate collection center at the end of its useful life. Selective and appropriate collection for the recycling of no longer used appliances and their disposal and destruction in an environmentally friendly manner, helps to prevent potential negative effects on the environment and health, and encourages the recycling of materials used to make the product.

For more detailed information about the treatment, recovery and recycling of this product, please contact your local authority, your household waste disposal service or the retailer where you purchased the product.

9. WARNING

As we continuously improve our products for the benefit of our customers, we reserve the right to modify technical characteristics without notice.

Warranties for MEDIFROID brand products are offered exclusively by selected retailers. No part of these instructions may be considered as a supplementary guarantee.

FRIO ENTREPRISE cannot be held liable for errors or technical or editorial omissions in this document.

Non-binding document.

Set-up conditions	
Degrees of pollution (Surrounding the appliance)	PD2
Ambient temperature	16~32°C
Maximum relative ambient moisture	80%
Maximum operating elevation (meters above sea level)	2000

Electrical values	PHAR200	PHAR300/400
Nominal voltage	220-240	220-240
Frequency	50	50/60
Connected load	0.98	2.0
Maximum mains fluctuation	±10%	±10%
Overvoltage categories	II	II

Lighting	
Energy efficiency class	Light source
This product contains one or more energy efficiency class F light sources	LED

Data according to DIN 13277				
Model	Energy consumption (kWh/24h)	Temperature increase-power failure (min.)	Temperature homogeneity (°C)	Temperature consistency (°C)
PHAR200	1,6	60	1,5	2,8
PHAR300	2,9	90	1,8	2,5
PHAR400	3,3	100	1,4	2,1

Appliance battery	PHAR200	PHAR300/400
Type	BT-12M1.3AT	BT-12M4AT
Nominal voltage	12V	12V
Capacity	1.3Ah	4.0Ah

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